

WARNING: Read and follow all directions carefully to avoid possible injury.

COLD USE Directions For Use:

1. Remove gel pack from the stuffed animal by unzipping the zipper on the side of the stuffed animal. 2. Store in the freezer for at least 1 hour prior to use. 3. Once frozen, remove gel pack from freezer, place gel pack in stuffed animal and place on affected area. Always use product with the provided cover. Never apply directly to skin. NOTE: This unit includes two (2) gel packs for convenience. Store one of them in the freezer for immediate access when needed.



HOT USE

BE EXTREMELY CAREFUL WHEN HANDLING PACK FROM MICROWAVE TO AVOID BURNS. HEAT ONLY IN A MICROWAVE. HEATING TIME MAY VARY DEPENDING ON MICROWAVE RATED OUTPUT.

HEATING INSTRUCTIONS

From Room Temperature: 1. Heat in microwave on HIGH in 30 second intervals until a safe and desired heat level is reached, with total heating time not to exceed 90 seconds. Massage and turn gel pack over in microwave between intervals to avoid hot spots. 2. Monitor gel pack in microwave while heating. If any swelling of the gel pack occurs, stop heating immediately and wait for the swelling to subside before removing. 3. Before removing gel pack from the microwave, carefully check for any signs of possible rupture or leakage. If any leakage is found, do not use the gel pack.

From Frozen: 1. Heat in microwave on HIGH for 60 seconds, then in intervals of 30 seconds until a safe and desired heat level is reached, with total heating time not to exceed 120 seconds. Massage and turn gel pack over in microwave between intervals to avoid hot spots. 2. Monitor gel pack in microwave while heating. If any swelling of the gel pack occurs, stop heating immediately and wait for the swelling to subside before removing. 3. Before removing gel pack from the microwave, carefully check for any signs of possible rupture or leakage. If any leakage is found, do not use the gel pack.

PLACE PRODUCT IN THE STUFFED ANIMAL COVER PRIOR TO USE.
NEVER APPLY DIRECTLY TO SKIN.

STORAGE: Keep gel pack in a freezer when not in use.

CAUTION 1. Consult a physician prior to use. 2. People with poor circulation, diabetes, nerve damage, paralysis or sensitive skin should not use this product. 3. Although non-toxic, discard gel pack if product leaks or is punctured. 4. Adult supervision is required for use on invalids or children. Do not use on infants or children under the age of 2 years. 5. Do not use the gel pack in conjunction with liniments, sports creams, balms or lotions. 6. Do not use on broken skin. 7. If gel pack is too hot, allow to cool before use. 8. Do not fall asleep with the gel pack.