IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY

DANGER: TO REDUCE RISKS OF BURNS, ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

1. READ ALL INSTRUCTIONS.
2. DO NOT USE WHILE SLEEPING.
3. DO NOT USE ON AN INFANT OR ON AN ANIMAL.
4. THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, A SLEEPING OR UNCONSCIOUS PERSON, A PERSON WITH DIABETES, OR A PERSON WITH POOR BLOOD CIRCULATION.
5. DO NOT USE ON AREAS OF INSENSITIVE SKIN.
6. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING. CHECK SKIN UNDER PAD FREQUENTLY.
7. **DO NOT USE IN AN OXYGEN ATMOSPHERE.**
8. **DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.**
9. **DO NOT SIT ON OR AGAINST OR CRUSH PAD-AVOID SHARP FOLDS.**
10. **NEVER PULL THIS PAD BY THE SUPPLY CORD.**
11. **DO NOT USE THE CORD AS A HANDLE.**
12. CAREFULLY EXAMINE INNER COVER BEFORE EACH USE. DISCARD THE PAD IF INNER COVERING SHOWS ANY SIGNS OF DETERIORATION.
13. NEVER USE PAD WITHOUT THE REMOVABLE COVER IN PLACE.
14. THIS UNIT SHOULD NOT BE USED BY CHILDREN WITHOUT ADULT SUPERVISION.
15. PLACE PAD ON TOP OF AND NOT UNDER THE PART OF THE BODY NEEDING HEAT.
16. USE THIS PAD ONLY ON A 110-120 VOLT AC CIRCUIT. UNPLUG WHEN NOT IN USE.
17. DO NOT TAMPER WITH THIS PAD IN ANY WAY.
18. THERE ARE NO USER-SERVICEABLE PARTS. IF FOR ANY REASON THIS PAD DOES NOT FUNCTION SATISFACTORILY, SEE WARRANTY FOR CONSUMER SERVICE INFORMATION.
19. DO NOT USE THIS PAD WITH LINIMENT, SALVE OR OINTMENT PREPARATIONS, ESPECIALLY ONES THAT CONTAIN HEAT-PRODUCING INGREDIENTS. SKIN BURNS COULD RESULT.
20. DO NOT BEND OR PINCH CORD AND DO NOT USE IF CORD IS DAMAGED.
21. **NEVER LEAVE THE APPLIANCE UNATTENDED, ESPECIALLY IF CHILDREN ARE PRESENT.**
22. **LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS.**
23. **SAVE THESE INSTRUCTIONS.**

This appliance has a polarized plug (one blade is wider than the other). This plug fits a polarized outlet only one way; if the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician.

DO NOT ATTEMPT TO DEFEAT THIS SAFETY FEATURE.

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FREQUENTLY ASKED QUESTIONS

1. **When is it best to use a heating pad?**
Use your heating pad whenever hot applications are desirable for personal comfort, and as recommended by your health care professional for heat therapy.
2. **How long does it take for the heating pad to warm up?**
This depends primarily on the heat setting you have selected. You will begin to feel the pad warming up within a minute or so.
By Setting: Warm – 3 Minutes Low – 6 Minutes Medium – 8 Minutes High – 10 Minutes
3. **What temperature does my heating pad reach?**
Maximum heat temperature will not exceed the Underwriters Laboratory temperature limits (176 F - 80 C) when used on the High setting. Temperatures vary by heat setting. You may want to use a lower heat setting depending on application and preference.
4. **Can I wash my heating pad cover?**
Yes, please see instructions inside this manual.
5. **Can I sit on or bend my heating pad?**
Do not sit on, bend or crush the heating pad. Avoid sharp fold as damage to wires could occur.

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IMPORTANT NOTES

Please read this entire manual carefully before using this device. Product may vary from picture.

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